

# Monitor Neurological Status During Any Conscious Procedure

*Choose the B-SMART®*

*Buchbinder Sensory Motor Activated Response Timer*

- *Monitor neurological function during conscious sedation for any procedure*
- *Monitor level of consciousness during administration of neuroleptic anesthesia*
- *Alerts OR staff when a patient is not responding to a vibratory stimulus*
- *Alerts OR staff to changes in a patient's neurological status simply, reliably, and accurately*
- *Timed repetitive signals can be set for 30, 60 or 120 second intervals*
- *Patient interactive cooperation is an added safety measure*
- *Easy set-up*



VASCULAR INSTRUMENTS  
**KOVEN**  
TECHNOLOGY INCORPORATED

**Koven Technology, Inc.**  
Your Vascular Healthcare Partner