

Diagnosis is the First Step

Choose the Bidop® 3 Vascular Doppler for early diagnosis of vascular disease

- Fully customizable - choose only the options & accessories to meet the specific needs of your facility
- Real-time waveform LCD - displays the waveform during testing
- 30 waveform memory storage - save waveforms to view later or download to Smart-V-Link® software
- Automatic gain and baseline control - no manual adjustments required
- Smart-V-Link® vascular software available for EMR's and documentation for reimbursement
- Flat vascular probe available - eliminates the need to hold the probe during testing
- PPG available to expand vascular testing capabilities
- Doppler rollcart available for user convenience
- Spanish display option



Smart-V-Link® vascular testing software is available for professional presentation & hard-copy reports for reimbursement without usage fees or significant hardware costs (run on any Windows® based computer).

*Convenient, handheld design ...
Results you can see*

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner