


*For Life & Limb*

**BSAVED**  
BY KOVEN

BUCHBINDER SURGEON ADJUSTABLE VALVE EXCISER DEVICE

- *Designed by a Vascular Surgeon to maximize the safety and effectiveness of valvulotomy procedures*
- *Enables the surgeon to maintain control by adjusting the blade diameter while inside the vessel*
- *Ergonomic handle design with single-handed blade release and retrieval*
- *Stationary blade deployment ensures that the cutting head stays in place as the blades are deployed*
- *Blades set into position for precise vein diameter adjustment*
- *Unique design prevents blades from coming in contact with the vein wall*
- *Blades capture valves while in closed position for more effective cutting*
- *Equally effective in small and large veins*
- *Pre-sterile, single-use*

**BSAVED Valvulotome**  
*For In-Situ Bypass Graft Procedures*



VASCULAR INSTRUMENTS  
**KOVEN**  
TECHNOLOGY CANADA INC.  
Koven Technology Canada Inc.  
Your Vascular Healthcare Partner