

Diagnosis is the First Step

Choose the Smartdop™ 45 Vascular Testing System for early diagnosis of vascular disease

- Convenient hand-held size includes features found in larger systems
- Fully customizable - choose only the options & accessories to meet the specific needs of your facility
- ABI and PVR studies for macrovessel assessment
- PPG Toe Pressures for microvessel assessment
- Unique, real-time waveform LCD - displays the waveform before printing
- 30 waveform memory storage - save waveforms to view later or download to Smart-V-Link software
- Automatic gain and baseline control - no manual adjustments are required
- Waveform strip generated for immediate documentation (or choose Smart-V-Link vascular software for full size professional reports & EMR's)
- PPG and PV modules available to expand testing capabilities
- French display option



*Unique,
handheld design...
Proven results*



Smart-V-Link vascular testing software is available for professional presentation & reporting without usage fees or significant hardware costs (run on any Windows based computer).

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY CANADA INC.

Koven Technology Canada Inc.
Your Vascular Healthcare Partner