

Diagnosis is the First Step

Choose the Smartdop™ 50EX-F bidirectional Doppler with frequency analysis for early diagnosis of vascular disease

- Frequency analysis mode included for carotid & penile studies
- ABI and PVR studies for macrovessel assessment
- PPG Toe Pressures for microvessel assessment
- Unique, real-time waveform LCD display - displays the waveform before printing
- 32 waveform memory storage - save waveforms to view later
- Automatic gain and baseline control - no manual adjustments are required
- Waveform strip generated for immediate documentation (or choose Smart-V-Link® vascular software for full size professional reports & EMR's)
- PPG and PV modules available for expanded testing capabilities



Smart-V-Link® vascular testing software is available for professional presentation & reporting without significant hardware costs (run on any Windows based computer) or usage fees

Real-time Doppler with frequency analysis capabilities

VASCULAR INSTRUMENTS
KOVEN
 TECHNOLOGY CANADA INC.

Koven Technology Canada Inc.
 Your Vascular Healthcare Partner