

Sensitive to your Surgical Probe Demands

- Durable and Cost-effective
- Surgical probes are interchangeable with Koven Surgical Dopplers
- Specialty probes for Vascular, Cardiovascular, Laparoscopic, Neurosurgery, Plastic and Reconstructive
- Autoclavable probes available
- Wide variety of sterilization options available (ETO, Steris, Sterrad, Autoclave)
- Non-invasive probes also available (2, 5, 8, and 10 MHz)
- 1 Year warranty *

* 1 year parts & labor warranty when used as recommended. Warranties for ACP, NRP, and CRP series probes differ. See warranty information for details.



Surgical Pencil Probes

Available in 5, 8, & 10 MHz
Sterilization: ETO, Sterrad



Autoclavable Probe

Model: ACP-8 (8 MHz)
Sterilization: Autoclave
* Adaptor required



Surgical Thin Probes

Models: VRP-8 & VRP-20 (8 & 20 MHz)
Sterilization: ETO, Sterrad, Steris
* Adaptor required



Cardiovascular Probes

Models: CRP-10H & CRP-20H (10 & 20 MHz)
Sterilization: ETO, Sterrad, Steris
* Adaptor required



Bayonet Style Neurovascular Probe

Models: NRP-10H & NRP-20H (10 & 20 MHz)
Sterilization: ETO, Sterrad, Steris
* Adaptor required



Laparoscopic Probe

Model: LRP-8 (8 MHz)
Sterilization: ETO, Sterrad, Steris
* Adaptor required



Flat Monitoring Probe

Model: FDP-8 (8 MHz)
Sterilization: ETO, Sterrad
* Adaptor required



Air Emboli Flat Monitoring Probe

Model: BF2M20S8A (2 MHz)
Compatible with ES-100VX and Bidop® 3 Noninvasive Dopplers



Koven Technology Canada Inc.
Your Vascular Healthcare Partner

T: 204.943.2596 / 204.943.2624 • F: 204.943.2652 • info@koven.com • www.koven.ca • email PO's: canadapurchase@koven.com