


# Monitor your patients... Even when you're not around

## *Choose the DS-250 Ambulatory Blood Pressure Monitor*

- Review patient response to blood pressure medication at work, home, and during sleep
- Verify time and duration of hypertensive episodes
- Dose titration of blood pressure medication
- Wide selection of time settings (12, 24 & 48 hour monitoring)
- Dual simultaneous measurements (Oscillometric [OSC] or Karotkoff [K-Sounds])
- Easy to use MS Windows® based software
- Flexible interval programming
- Compact monitor & cuffs for portability and user convenience
- USB connection for convenient uploading & downloading



### Data Samples



Trend



Correlation

VASCULAR INSTRUMENTS  
**KOVEN**  
TECHNOLOGY INCORPORATED

**Koven Technology, Inc.**  
Your Vascular Healthcare Partner