

Sensitive to your Surgical Doppler Demands

Choose the ES-100X Surgical MiniDop®

- **Wide selection of interchangeable probe styles and frequencies for user flexibility**
- **Autoclavable and single-use probes available**
- **Specialty probes for Vascular, Laparoscopic, Plastic & Reconstructive, Cardiovascular, and Neurosurgery**
- **Surgical and non-invasive probes available in 5, 8, 10 and 20 MHz**
(non-invasive available in 5, 8 & 10 MHz)
- **1-year warranty on Doppler and probes reduces replacement costs**
(See warranty information for CRP, NRP, and ACP series probes)
- **Compact and battery operated**
- **Doppler rollcart and footpedal remote available for convenience & surgeon control in the OR**



Footpedal Remote

Doppler Rollcart

Optional accessories include a Doppler rollcart and footpedal remote

*Compact & Convenient...
Ideal for the OR*




Koven Technology, Inc.
Your Vascular Healthcare Partner



Pencil Probe
(P8M0558A 8 MHz)
* also available in 5 & 10 MHz



Single-Use
(VDP-8 8 MHz)



Surgical Specialty Probes
(8, 10, & 20 MHz available)
Cardiovascular, Laparoscopic, Neurovascular
and autoclavable models available

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • info@koven.com • www.koven.com • email PO's: purchase@koven.com