


Sensitive to your Surgical Doppler Demands

Choose the HD-307 Bidirectional Surgical Doppler

- **Wide selection of interchangeable probe styles for user flexibility**
- **Autoclavable and single-use probes available**
- **Specialty probes for Vascular, Laparoscopic, Plastic & Reconstructive, Cardiovascular, and Neurosurgery**
- **Surgical probes available in 8, 10 and 20 MHz**
- **1-year warranty on Doppler and probes reduces replacement costs**
(See warranty information for CRP, NRP, and ACP series probes)
- **Single-Frequency Bidirectional Doppler**
(8, 10 & 20 MHz models available)
- **Compact and battery operated**
- **Doppler rollcart and footpedal remote available for convenience & surgeon control in the OR**




*Lightweight & Portable ...
Bidirectional Flow Assessment*



Pencil Probe
(PP8M05N5A 8 MHz)
* also available in 10 MHz



Single-Use
(VDP-8 8 MHz)



Surgical Specialty Probes
(8, 10, & 20 MHz available)
Cardiovascular, Laparoscopic, Neurovascular
and autoclavable models available

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner