


Diagnosis is the First Step

Choose the *Smartdop® 30EX Vascular Testing System*
for early diagnosis of vascular disease



*Durable & Effortless...
Automated, reliable
vascular testing*

- **Flat vascular probe eliminates the need to hold the probe during testing**
(pencil style probes also available)
- **Fully customizable - choose only the options & accessories to meet the specific needs of your facility**
- **Built-in automatic cuff inflator - quickly performs automated ABI, PVR and PPG toe pressure studies**
- **Automatically calculates and prints ABI and TBI results**
- **Unique, real-time waveform LCD - displays the waveform before printing**
- **30 waveform memory storage - save waveforms to view later or download to Smart-V-Link® software**
- **Automatic gain and baseline control - no manual adjustments required**
- **Waveform strip generated for immediate documentation for reimbursement**
- **Smart-V-Link® vascular software available for EMRs**
- **Optional PPG probe available for expanded testing capabilities**
- **Doppler rollcart available for user convenience**
- **Spanish display option**



Smart-V-Link® vascular testing software is available for professional presentation & reporting without usage fees or significant hardware costs (run on any Windows® based computer).

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner