

# *Sensitive to your Surgical Probe Demands*

- Durable and Cost-effective
- Surgical probes are interchangeable with Koven Surgical Dopplers
- Specialty probes for Vascular, Cardiovascular, Laparoscopic, Neurosurgery, Plastic and Reconstructive
- Pre-sterile single-use & autoclavable probes available
- Wide variety of sterilization options available (ETO, Steris®, Sterrad®, Autoclave)
- Non-invasive probes also available (2, 5, 8, and 10 MHz)
- 1 Year warranty \*

\* 1 year parts & labor warranty when used as recommended. Excludes VDP-8 single-use. Warranties for ACP, NRP, and CRP series probes differ. See warranty information for details.

## **Cardiovascular Probes**

Models: CRP-10H & CRP-20H (10 & 20 MHz)

Sterilization: ETO, Sterrad®, Steris®

\*Adaptor required

## **Bayonet Style Neurovascular Probe**

Models: NRP-10H & NRP-20H (10 & 20 MHz)

Also available in **Flexible Tip** models HP-10 & HP-20 (10 & 20 MHz) for easier arterial access to difficult locations

Sterilization: ETO, Sterrad®, Steris®

\*Adaptor required

## **Laparoscopic Probe**

Model: LRP-8 (8 MHz)

Sterilization: ETO, Sterrad®, Steris®

\*Adaptor required

## **Flat Monitoring Probe**

Model: FDP-8 (8 MHz)

Sterilization: ETO, Sterrad®

\*Adaptor required

## **Air Emboli Flat Monitoring Probe**

Model: BF2M20S8A (2 MHz)

Compatible with ES-100VX and Bidop® 3 Noninvasive Dopplers

## **Surgical Pencil Probes**

Models: P5M05S8A, P8M05S8A, P10M05S8A (5, 8 & 10 MHz)

Models for HD-307: PP8M05N5A, PP10M05N5A (8 & 10 MHz)

Sterilization: ETO, Sterrad®

## **Autoclavable Probe**

Model: ACP-8 (8 MHz)

Sterilization: Autoclave

\* Adaptor required

## **Surgical Thin Probes**

Models: VRP-8 & VRP-20 (8 & 20 MHz)

Sterilization: ETO, Sterrad®, Steris®

\*Adaptor required

## **Single-Use, Pre-Sterile Probes**

Model: VDP-8 (8 MHz)

\*Adaptor required



**Koven Technology, Inc.**  
Your Vascular Healthcare Partner

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • [info@koven.com](mailto:info@koven.com) • [www.koven.com](http://www.koven.com) • email PO's: [purchase@koven.com](mailto:purchase@koven.com)