


For Life & Limb

The Koven Valvulotome

For In-Situ Bypass Graft Procedures

- **Designed by a Vascular Surgeon to maximize the safety and effectiveness of valvulotomy procedures**
- **Enables the surgeon to maintain control by adjusting the blade diameter while inside the vessel**
- **Ergonomic handle design with single-handed blade release and retrieval**
- **Stationary blade deployment ensures that the cutting head stays in place as the blades are deployed**
- **Blades lock into position for precise vein diameter adjustment**
- **Unique design prevents blades from coming in contact with the vein wall**
- **Blades capture valves while in closed position for more effective cutting**
- **Equally effective in small and large veins**
- **Pre-sterile, single-use**



VASCULAR INSTRUMENTS
KOVEN[®]
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner